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1. Hexagonal screwdriver, size 3.0 mm - general use, mainly various flange screws, the most common size used for screws.
2. Flat screwdriver, size 3.2 mm - general use.
3. Vacuum grease - used for sealing of the o-rings.
4. The basic specimen stubs.
5. Carbon conductive adhesive discs for bonding the sample to a stub.
6. Hexagonal screwdriver, size 2.5 mm - used for mounting / dismantling sample stage Z extensions (LS, LM chambers) and the stub holder (XM chamber).
7. Hexagonal screwdriver, size 1.5 mm - used for fixing the screws in the sample stage.
8. Flat tweezers - general use.
9. Aperture 75 μ m.
10. Low vacuum aperture holder.
11. Screw, size M3x25 - used for the removal of the low vacuum aperture holder.

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
The **standard specimen holder diameter 12.5 mm, height 3 mm** - The most common holder suitable for any smaller specimens. The specimen is normally glued by glue or sticked on the double side sealing tape.


The **specimen holder HM-114-R** - The holder is suitable for flat samples up to the width of 16 mm.

The **specimen holder HM-116-R** - The holder is suitable for flat samples up to the width of 40 mm.


The **specimen holder ACB-16-00** - The holder is suitable for flat samples up to the width of 4 mm.


The **specimen holder GM-163-R** - diameter 12 mm, height 12 mm. The extension for small specimen holders to be able to work on very short working distances.

The **specimen holder TE00038R** - The holder is designed for standard samples of 30 mm diameter.

The **specimen holder GM-153-R** - The replacement of standard seven position specimen holder (the seven position holder must be removed) it is designed for rod shaped samples of diameter up to 26 mm.

The **specimen holder GM-154-R** - CAMSCAN 12.5 mm specimen holder adapter.

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Info Panel		
Continual	Single	Acquire
Scan Mode:	RESOLUTION	
HV:	30.00 kV	
Magnification:	40.00 kx	
View field:	3.611 μm	
Speed:	1 (0.10 μs/pxl)	
WD:	5.000 mm	
Depth of Focus:	1.07 μm	
Stigmator:	0.1 % / 0.2 %	
Shift:	0.0 μm / -0.0 μm	
Rotation:	0.00 deg	
Beam Intensity:	9.00	
Absorb. Curr:	180 pA	
Spot Size:	2.4 nm	

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
Function name	Displayed value
Preset	Name of the actual preset
Beam energy	The actual energy of the ion beam
Condenser	Voltage on condenser electrodes
Probe aperture	The index and diameter of the current probe aperture diameter of the ion beam.
Magnification	The actual magnification of the image (the ratio of the scanning window width to the real width of the scanning area on the sample surface).
Speed	The index of the current scanning speed, the value in brackets shows current time per pixel
Image shift	The actual values of the electronic image shift
Rotation	The actual values of electronic rotation of the image
Stigmators	The actual astigmatic correction values
Probe current	The beam current measured in the bottom part of the column - after the beam has passed the condenser and apertures. The value is valid only if the ion beam is blanked (the scanning is stopped). The measured current is equivalent to the real beam current impinging on the specimen during ion beam scanning.
Objective	Displays the current voltage on the objective electrodes
Gain Blacklevel	/ Displays the actual values of the detector gain and offset

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
Object Panel		Object Edit Panel	
	Select and edit objects		Show/hide background image. To do so, select and click (left mouse button) into image window of your interest.
	Create a point object		Show/hide grid (If grid is enabled, the snapping is active.)
	Create a line object		Create a drift correction mark
	Create a cross object		Show/hide all objects in the project
	Create an outline rectangle		Undo
	Create a filled rectangle		Redo
	Create a polishing rectangle		Cut
	Create a stairs rectangle		Copy selected objects
	Create a reference point object		Paste selected objects
	Load a bitmap or another format object		Delete selected objects
	Insert text		Clone selected objects
	Create an outline polygon		Align selected objects
	Create a filled polygon		Group selected objects (Objects can be grouped only within one layer.)
	Create an outline circle		Break a group of objects
	Create a filled circle		Create mesh of selected object
	Create a circle annulus		Rotate selected object counter-clockwise by 90 degrees
	Create a polishing circle		Rotate selected object clockwise by 90 degrees
	Create a stairs circle		Create a void rectangle
	Create a void circle		Create a void polygon

Gas Injection System界面功能介绍


Heating or Cooling - click this button to start heating or cooling of the reservoir, capillary and the nozzle (depends on the precursor type). When the temperatures reach the values required for operation, the bar graphs above become green and the Valve button is enabled. Click this button again to finish using the precursor and turn the heating off.

Outgas - click this button to start the outgas procedure manually. This procedure must be performed, when the reservoir has not been used for long period of time or when the system has been vented. If the reservoir is properly degassed, the outgas procedure should not increase the pressure in the chamber - in such a case we can interrupt the outgas procedure by pressing the button again.

Note: For some precursors e.g. W, the outgas procedure is forced during turning off the heating - to eliminate nozzle clogging.



Valve - click this button to open and close the precursor's valve. It is only allowed to open the valve when the temperatures of the reservoir, capillary and the nozzle are on the values required for operation. By opening the valve, you start the actual flow of the gas through the capillary and nozzle to the microscope chamber.

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Nanomanipulator界面功能介绍

For the most typical 3D nanomanipulator, there are three axis of movement (XYZ). The orientation of the axis is set to be the same as the **specimen stage** i.e. +X moves to the right, +Y moves up, +Z moves down.

- By holding the **arrow** button, manipulator moves with the selected speed
- **X,Y,Z** shows the current manipulator position in absolute coordinates in μm (position from the calibration mark)
- Speed defines the speed of the movement for the manual movement. SCAN1 and SCAN2 speed indicates the mode where the manipulator moves by changing the length of piezo elements only, with no jumping and vibrations, but with movement range limited to $\sim 1 \mu\text{m}$ only. When entering the absolute coordinate or selecting the saved position it will move at the maximum speed.
- The **Calibrate** button starts the automated calibration procedure.
- There are three basic **Positions** selectable.
 - **Parking position** is set to be a safe position for any stage movement without a risk of collision. This position is set during the installation and can not be changed by the user. In this position also the position sensors are switched off, so there should be no interference to EDX, EBSD, CL or other detectors.
 - **Working position** is defined by the user and **is stored for each user separately**.
 - **Standby Position** is defined relatively to the working position (by default $200 \mu\text{m}$ above the working position).

The movement from parking to working position (or back) is always done through the Standby position.

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